Smart Implant Solutions

SIS DIGITAL ANALOG PROCEDURE

1. SIS digital analog calibration
2. Step by Step

PROCEDURE

1. SIS digital analog

PREVIOUS:

In each **Modelcreator library** you will find a folder called **“Patrones fricciones”**.

Within this folder you will have **several folders depending on the sizes** of each Implant system.

Example:

1 Copy digital analog libraries to path:

GALWAY 3.0 > exocad-DentalCAD3.0-2021-03-25 > DentalCADApp > library > modelcreator > implants

The image displays three sequential screenshots of a Windows File Explorer window, illustrating the directory structure for digital analog libraries. The first screenshot shows the root directory 'implants' with a list of folders including 'SIS-0900_NB-AC'. The second screenshot shows the 'SIS-0900_NB-AC' folder containing subfolders 'fricciones', 'Analog', 'Analog2', 'Base', 'Base2', and 'config'. The third screenshot shows the 'fricciones' folder containing subfolders 'NP' and 'RP'. Arrows indicate the navigation path from the root to the specific folder and then to the subfolders.

2 Within each folder **patterns** and **frictions** are displayed.

3 Within the **“Fricciones”** folder, the available sizes are displayed.

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
4 The "fricciones" folder contains:

- 5 subfolders with the 5 frictions.
- The pattern of your system.

5 Each subfolder contains:

- The file with each friction.

... implants > SIS-0900_NB-AC > Fricciones > NP > 3

PROCEDURE

2. Step by Step

1 Select the print pattern.

2 Print with the resin for models at the usual resolution.

4 Try the analog in all frictions and **choose** the friction.

3 Wash and **Cure** as usual.


PROCEDURE


2. Step by Step

5 Go to the **subfolder** and **copy** the file.

... implants > SIS-0900_NB-AC > Fricciones > NP


... implants > SIS-0900_NB-AC > Fricciones > NP > 3

Copy file

6 Finally, go to the **root folder** of the library and **paste the file**.

... implants > SIS-0900_NB-AC

Paste file

The file will be replaced automatically.

Your library is ready to be used with the chosen friction.
Don't forget to **restart the software**.

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