DESCRIPTION

Screw destined to fix any prosthetic element to the implant.

HOW TO USE / RECOMMENDATIONS

It is recommended not to use screws from the prosthetic laboratory in clinical use. For its screwing/unscrewing it is important to use the appropriate size wrench. In case of doubt, check that the upper measuring key does not fit into the housing.

The torque recommended by Smart Implant Solutions for permanent prostheses, unless otherwise stated on the product label, is as follows:

METRIC	TORQUE	METRIC	TORQUE
M2,5	30-35 N*cm	M1,6	25-30 N*cm
M2	30-35 N*cm	M1,4	15-20 N*cm
M1,8	25-30 N*cm		

MATERIAL

Grade 5 titanium.

FINISH


Optionally, they can be supplied with **BIOCARB** coating.

ADVANTAGES

They stand out for their great hardness and mechanical resistance.


BIOCARB:

- Amorphous carbon multilayer coating.
- High tenacity and hardness (2000-5000HV).
- Gray finish.
- Fully heme and biocompatible as it is carbon based.
- Sterilizable.
- High resistance to wear.
- Better screw tightening.
- Reduction of the torque applied to the screw.
- The implant screws can be removed without risk of breakage or seizure.

DATA SHEET

LAB SCREW SMART ANGLE


DESCRIPTION

Screws with Smart Angle connection that allows screwing at an angle from 0° to 30°. Once the crown or bridge has been fabricated in the laboratory, it is fixed directly to the implant using the Smart Angle screw.

HOW TO USE / RECOMMENDATIONS

It is screwed on an interface and on a Cr-Co base, using the Smart Angle wrench. *For the recommended tightening torques, see the technical data sheet for screws.*

SMART ANGLE SCREWDRIVER TIP

SHORT	MEDIUM	LONG
20mm	25mm	35mm
SIS-9142C	SIS-9142M	SIS-9142L

MATERIAL

Grade 5 titanium.

FINISH

Blue anodized.

ADVANTAGES

- Safest angled solution on the market.
- Its improved connection avoids crowning of the screws.
- The interfaces, Cr-Co bases and implant libraries are ready to use the Smart Angle angled solution.

Blue anodizing allows distinguishing a standard screw from a Smart Angle screw.

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